PHILIPS

Transcription headphones

Under-the-chin style

Gold-plated plug

14 mm speaker driver

3 meter/10-foot cable

LFH0334



Outstanding audio quality

with Philips transcription headphones

Lightweight under-the-chin style stereo headphones designed to deliver excellent sound quality, with soft ear cushions for wearing comfort, and a hanging bracket that attaches the headset to a monitor.

Superior audio quality

- A 14 mm speaker driver optimizes high audio quality
- Gold-plated plug ensures an ultra reliable connection
- Neodymium magnet enhances sensitivity and sound quality

For optimal results

- + Lightweight design improves comfort for prolonged use
- Under chin style to fit all ears for better comfort
- \cdot 3 meter/10-foot cable for freedom to move about your work area



14 mm/half-inch speaker driver

Small enough for optimum wearing comfort, yet large enough to deliver crisp, nondistorted sound, the 14 mm/half-inch speaker driver is the ideal size to deliver high quality sound.

Gold-plated plug

The precious-metal gold finishing on the plug will produce a more reliable connection for better audio quality.

Neodymium magnet

Neodymium is the best material for producing a strong magnetic field for greater sensitivity and higher overall sound quality.

Lightweight design

The durable, lightweight quality materials enhance comfort for extended wearing.

Under chin style

Designed for an under-the-chin fit for better comfort.

3 meter/10-foot cable

The 3 meter/10-foot cable gives you the freedom to move around your work area without having to remove the headphones.

Specifications

Sound

Acoustic system: open Magnet type: neodymium Diaphragm: PET dome Voice coil: copper Sound quality: optimized for voice reproduction Impedance: 2 × 32 Ohm Maximum power input: 2 × 10 mW Sensitivity: 103 dB Speaker diameter: 14 mm Type: dynamic System: stereo

Connectivity

Cable connection: one-sided Cable length: 3 m/10 ft Connector: 3.5 mm stereo, right angled Connector finish: gold-plated Type of cable: copper

Package contents

Headphones Hanging bracket

Compatibility

All Philips desktop dictation and transcription systems Philips Pocket Memo digital voice recorder Philips Pocket Memo meeting recorder

Green specifications

Compliant to 2011/65/EU (RoHS) Lead-free soldered product

Specifications

Product dimensions: 160 mm × 150 mm × 18 mm/6.3" × 6.0" × 0.8" Weight: 53 g/1.9 oz

Accessories Philips USB audio adapter LFH9034

© 2021 Speech Processing Solutions GmbH. All rights reserved.

Specifications are subject to change without notice. Philips and the Philips shield emblem are registered trademarks of Koninklijke Philips N.V. and are used by Speech Processing Solutions GmbH under license from Koninklijke Philips N.V. All trademarks are the property of their respective owners.



Issue date 2021-01-04 Version 1.0

www.philips.com/dictation